



## Product Data Sheet ZP50-240-9 Serial # 1263.3-1

Thank you for your purchase of an Xradia zone plate optic. In this document you find the parameters and fabrication tolerances of the zone plate.

### Zone Plate Parameters:

Outer Diameter	μm	240
Inner Diameter (no zones)	μm	N/A
Outermost Zone Width	nm	50
Zone Material		Electroplated Gold
Zone Height	nm	900 +/- 8%
Number Of Zones <sup>2</sup>		1200
Suggested Energy Range <sup>3</sup>	keV	5-24
Theoretical Max. Diffraction Efficiency <sup>1</sup>	%	17 (@ 9 keV)
Support Membrane Material		Si <sub>3</sub> N <sub>4</sub>
Support Membrane Thickness	μm	0.3
Support Membrane Size	mm	0.5 x 0.5
Support Silicon Frame Size	mm	2 x 2 x 0.18
Central Stop Diameter <sup>4</sup>	μm	80
Central Stop Height / Material	μm	85 / Au

<sup>1</sup> Xradia delivers zone plates with a focusing efficiency exceeding 50% of the theoretically calculated value.

<sup>2</sup> Number of fabricated rings calculated for a full zone plate (no missing inner zones).


<sup>3</sup> Energy range for which the theoretical efficiency is greater 10%.

<sup>4</sup> The placement accuracy of the stop is 10μm (typical 5 μm).

### Zone Plate Focal Length

The focal length of the zone plate is given by:

$$\begin{aligned}
 f &= \frac{OD \Delta R_n}{\lambda} \\
 &= \frac{12.0 \text{ mm}}{\lambda \text{ [nm]}} \\
 &= 9.678 \text{ mm gE [keV]}
 \end{aligned}$$



4075 Sprig Drive

Concord CA

94520-8535

tel 925 288-1818

fax 925 288-0310

www.xradia.com



## Product Data Sheet ZP50-240-9 Serial # 1263.3-1


### ***Use In Vacuum***


Zone plates can be used in vacuum. Several precautions need to be taken:

- No pressure differential can be applied between both sides of the zone plate wafer.
- Change the ambient pressure slowly. Avoid shocks generated by opening large conductance pneumatic valves.
- Zone plates can be damaged by ultra-sonic particles generated during venting of the vacuum system. Vent slowly and with a clean gas (use filter if possible). Avoid venting through a port which is in direct line-of-sight of the zone plate. Install a shield if necessary.
- Make sure that the zone plate heat load is small and removed through proper thermal connections. If in doubt, consult Xradia for advice.

### ***SEM Micrograph***

Outermost zones of the zone plate





Zone Plate Mounting Diagram

---

4075 Sprig Drive

Concord CA

94520-8535

---

tel 925 288-1818

fax 925 288-0310

[www.xradia.com](http://www.xradia.com)